TRABALHOS LABORATORIAIS / LABORATORIAL LAB SESSIONS

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HYDROLYSIS OF ETHYL ACETATE

Kinetic study and Reactor Dimensioning

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Introduction

The saponification reaction of ethyl acetate ($\text{CH}_3\text{COOCH}_2\text{CH}_3$, EtOAc) with sodium hydroxide (NaOH) can be described by the following stoichiometric equation:

The analytical method used to follow the evolution of the concentration of sodium hydroxide is the direct measure of the pH of the reaction mixture.

The reaction will be studied in the three lab sessions, with a growing degree of complexity. In the 1st lab session, the kinetic law will be studied using the initial rate method, estimating the global order of the reaction and the rate constant at room temperature. In the 2nd lab session, a new variable will be introduced: the temperature. The values of the rate constant and equilibrium constant will be estimated for two temperatures allowing to determine the activation energy, E_a , and the heat of reaction, ΔH ; for this determination, data from the 1st lab session will be used so that there is a set of data from three temperatures for the calculation of E_a and ΔH , respectively. In the 1st and 2nd lab sessions, a batch reactor will be used, in which each trial occurs under isothermal conditions. In the 3rd lab session (occurring simultaneously with the 2nd one), a plug flow reactor, PFR, is used, operating under non-isothermal and non-adiabatic conditions. The values of the rate constant, equilibrium constant, activation energy and heat of reaction, estimated with data acquired in the 1st and 2nd lab sessions, will be used to determine the heat transfer coefficient of the material that constitutes the PFR.

The data processing it will also be of increasing degree of complexity: in the 1st lab session, only molar balance equations for the batch reactor and kinetic law for irreversible reaction will be used; in the 2nd lab session will be used the molar balance equations of the batch reactor, the kinetic law for a reversible reaction, Arrhenius law and van't Hoff's law; finally, in the 3rd lab session, will be used the molar balance equations for the PFR reactor, kinetic law for the reversible reaction, Arrhenius law, van't Hoff's law and energy balance equation. As the theoretical expressions of the variation of the degree of conversion (2nd and 3rd lab sessions) and temperature (3rd lab session) are differential equations, its numerical resolution implies to use an approximation method. The derivative (continuous domain) will be replaced by its algebraic approximation (discrete domain). For this purpose, the Euler method will be used, which consists in finite difference approximations; the smaller the interval considered in the discretization, the more precise the approximation will be.

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DETERMINATION OF THE GLOBAL REACTION ORDER BY THE INITIAL RATE METHOD

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Abstract and objectives

The purpose of these lab sessions is the kinetic study of the system ethyl acetate/sodium hydroxide, by monitoring the time evolution of the NaOH concentration (through pH measurements). The global reaction order and the rate constant will be determined applying the initial rate method, taking from literature the partial orders in relation to each reagent. An isothermal batch reactor is used. For the data treatment the molar balance equation to the batch reactor will be combined with the kinetic law assuming the irreversibility of the reaction in the initial times.

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Attention: Please check if the pH meter resolution defined as 0.001 (3 decimal digits)

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Add to the reactor **160 mL** of water and **20 mL** of NaOH (0.05 M). Switch on stirring and allow pH to stabilize.

Start reaction by adding **20 mL** de AcOEt (0.05 M) to the reactor and start counting time taking

t=0 at half addition.

Read pH every **10 s** during **7 min** and afterwards, read every **1 min**, up to **15 min**. Use in the pH meter the option: measuring continuously.

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Add to the reactor **130 mL** of **water** and **35 mL** of **NaOH** (0.05 M). After pH stabilization, add **35 mL** of **AcOEt** (0.05 M) and start counting time and pH measurements as before (at half addition).

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Add to the reactor **100 mL** of **water** and **50 mL** of **NaOH** (0.05 M). After pH stabilization add **50 mL** of **AcOEt** (0.05 M) and start counting time and do the pH measurements as before.

NOTAS/NOTES:

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Pay attention to the following aspects:

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- As the kinetic law used for data processing assumes an irreversible reaction, which occurs in the initial times, it is important that the pH meter provides an immediate response to a change in the concentration of the NaOH reagent. In order to figure out the response time of the pH electrode you are using, measure the time that elapses between the NaOH addition and the stabilization of the pH reading.

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- Determine which pH value the electrode should measure after adding the NaOH volume at the beginning of each trial (3 calculations).

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- During each test, take note of room temperature.

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- The same pH meter should be used in both lab sessions, 1 and 2, therefore, note the number of the electrode that is on your laboratory bench in lab session 1, so that you can use it again in lab session 2.

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[1] - Petek, A. and Krajnc, M., The enthalpy and entropy of activation for ethyl acetate saponification. Int. J. Chem. Kinet., **2012**, 44: 692-698. <https://doi.org/10.1002/kin.20712>

[2] - Das, K., Sahoo, P., Sai Baba, M., Murali, N., & Swaminathan, P. Kinetic studies on saponification of ethyl acetate using an innovative conductivity-monitoring instrument with a pulsating sensor. International Journal of Chemical Kinetics, **2011**, 43(11), 648-656

[3] - J. Murdoch, J.A. Barnes, Statistical Tables: for Students of Science Engineering Psychology Business Management Finance, Bloomsbury Publishing, **1998**, ISBN – 9780333558591.

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DETERMINATION OF ACTIVATION ENERGY AND HEAT OF REACTION

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Abstract and objectives

The purpose of this work is to obtain an estimation of the activation energy and heat of reaction of the saponification process, thus, there will be carried out experiments at different temperatures. The degree of conversion at each temperature will be simulated by the equation obtained from the condensation of the equation of the mass balance of the reactor with the kinetic law for a reversible second order reaction. The adjustment will be made by variation of the parameters of the kinetic constant, k_d , and equilibrium constant, K_e , having the minimization of the square of the difference of the simulated and the experimental values as criteria. As the theoretical expression of degree of conversion is a differential equation, the Euler method will be used to numerically estimate it. This method consists of approximation by finite differences, the approximation being more precise the smaller the time interval, Δt (domain of the variable monitored in this work), considered in the discretization.

The activation energy, E_a , will be estimated from the variation of k_d with the inverse temperature, according to the Arrhenius law:

$$k = k_{T_{Ref}} e^{\left(-\frac{E_a}{R} \left(\frac{1}{T} - \frac{1}{T_{Ref}}\right)\right)} \quad (\text{eq. 2})$$

The heat of the reaction, ΔH , will be estimated from the variation of K_e with the inverse temperature, according to Van't Hoff law:

$$K_e = K_{eT_{Ref}} e^{\left(-\frac{\Delta H_r}{R} \left(\frac{1}{T} - \frac{1}{T_{Ref}}\right)\right)} \quad (\text{eq. 3})$$

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Attention: *Please pay attention if the pH meter resolution defined as 0.001 (3 decimal digits)*

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Add to the reactor **197.6 mL** of **water** and **1 mL** of **NaOH** (1.0 M). Switch on the thermostatic bath and allow temperature to stabilize. When the reaction mixture reaches the desired temperature and the pH reading has stabilized, start reaction by adding **1.4 mL** of **AcOEt** (0.7 M) to the reactor. Start measuring the pH of the reaction mixture at the time corresponding to the half addition of ethyl acetate (t=0) and read the pH of the reaction mixture **every 10 s for 10 min**. Then read every **30 s until 20 min**.

Use in the pH meter the option: measuring continuously.

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In each of the tests, write down the temperature each time you record pH values. For the data processing class, organize your worksheet according to the following instructions.

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Lab session 2 - Preparation of the worksheet for the treatment of data

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$$X = \left(\frac{C_{Ao} - C_A}{C_{Ao}} \right) \quad (\text{eq. 4})$$

pH measurements

Register the time and pH values in a Excel lab worksheet. For data treatment, organize the lab worksheet according to Table 2, calculating pOH, $[\text{OH}^-]$ and conversion, X.

$$X = \left(\frac{C_{Ao} - C_A}{C_{Ao}} \right) \quad (\text{eq. 4})$$

Tabela 2 – Exemplo de registo inicial de tempo e medidas de pH; cálculo de pOH, $[\text{OH}^-]$ e X.

Table 2 – Example of initial register of time and pH values; calculation of pOH, $[\text{OH}^-]$ e X.

time	pH	pOH	OH^-	X_{exp}
0	11.455	2.36181	0.004347	0
10	11.45	2.36681	0.004297	0.011447
20	11.449	2.36781	0.004287	0.013721
30	11.442	2.37481	0.004219	0.02949
40	11.44	2.37681	0.004199	0.033949
50	11.439	2.37781	0.00419	0.036171
60	11.434	2.38281	0.004142	0.047204
70	11.426	2.39081	0.004066	0.064594
80	11.424	2.39281	0.004048	0.068892
90	11.42	2.39681	0.00401	0.077429
100	11.414	2.40281	0.003955	0.090087
110	11.413	2.40381	0.003946	0.092179
120	11.407	2.40981	0.003892	0.104635

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[1] - Petek, A. and Krajnc, M., The enthalpy and entropy of activation for ethyl acetate saponification. Int. J. Chem. Kinet., **2012**, 44: 692-698. <https://doi.org/10.1002/kin.20712>

[2] - Das, K., Sahoo, P., Sai Baba, M., Murali, N., & Swaminathan, P. Kinetic studies on saponification of ethyl acetate using an innovative conductivity-monitoring instrument with a pulsating sensor. International Journal of Chemical Kinetics, **2011**, 43(11), 648-656. <https://doi.org/10.1002/kin.20597>

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DETERMINATION OF THE HEAT TRANSFER COEFFICIENT OF THE WALL OF A NON-ISOTHERMIC AND NON-ADIABATIC PFR

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
Abstract and objectives:

Once the kinetic parameters, activation energy and heat of reaction, are determined in the previous works, one should obtain the heat transfer coefficient of the wall of a tubular reactor (PFR), working in non-isothermal and non-adiabatic conditions. The variation of temperature and degree of conversion will be determined, along the length of the reactor. The heat transfer coefficient, U , corresponds to the value that minimizes the difference between the values of the simulated and real temperatures. The equations that theoretically describe the temperature along the reactor and the degree of conversion are obtained by condensing the molar balance to the PFR reactor, with the kinetic law and the energy balance. As both theoretical expressions of degree of conversion and temperature variation are differential equations, the Euler method will be used to estimate them numerically. The approximation will be more precise the smaller the spatial interval, Δz , considered in the discretization. The temperature and conversion are monitored all along the reactor length, the z domain.

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Caption: (1) NaOH solution, (2) EtOAc solution, (3) pump, (4) heating bath, (5) PFR

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Volume do reator: A DETERMINAR

Reaction conditions:

- Concentrations of feed solutions: $C_{\text{NaOH}}=0.005 \text{ M}$, $C_{\text{EtOAc}}=0.1 \text{ M}$
- Feed volumetric flow rates:

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When the work is started, the heating bath should be already stabilized at **T=90°C** and the peristaltic pump working at a given **rate**. Write down that value (rpm).

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Start registering the **pH** and **reactor temperature** values **every minute** up to 10 min (ideally until steady state).

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Take note of **room temperature** and the **internal diameter** of the tube, on a piece of tube similar to that used in the set-up; take 5 measurements on each side and make an average. Measure the **length** of the tube in the set-up, between the electrodes.

Perform several **flow rate measurements** (ml/min), counting the amount of effluent collected in a graduated cylinder (ml) for 1 minute. Repeat 5 times and do an average.