
**Problem 1.** What would the expected effect be on a PCR reaction if the primers used were shorter than the intended oligonucleotide sequences?

- a) The PCR reaction would not commence
- **b)** The PCR reaction would end after one cycle
- c) The reaction would generate a single short PCR product
- d) The reaction would yield a mixture of non-specific products

### Problem 2.

Consider the DNA sequence below.

- a. Design primers (16 nts long) that will allow to obtain a PCR product with 400 pbs.
- b. Design a PCR program using the Taq enzyme (1000 pbs/min).



**Problem 3.** A researcher was trying to design a PCR experiment. The sequence he wished to amplify is part of a gene *murF* of *Staphylococcus aureus*, that is given below.

- 5' GCATGTCAATTGGGCCTATATGGCGTAGCAATTTGGCCGGCTATATGGCCGTAGC 3'
- 3' CCTACAGTTAACCCGGATATACCGCATCGTTAAACCGGCCGATATACCGGCATCG 5'

The researcher decided he could only afford to use primers of 6 bases in length.

- a) Write down the sequence of the two primers he should order to be made.
- b) The scientist carried out the experiment, using appropriate primers of 6 bases in length and a plasmid with the gene *murF* cloned. He then separated the pieces of DNA he had amplified. What technique would he use to separate the amplified pieces of DNA?
- c) The scientist obtained **two** bands of amplified DNA. Explain why at least two bands of amplified DNA would be obtained in the above experiment.
- d) What is the size of the 2 bands?
- e) When the scientist carried out the experiment using as template genomic DNA from *Staphylococcus aureus*, he counted 8 bands. Explain why he obtained so many more bands.
- f) Suggest how he could improve this experiment to obtain more reliable results.

# **Problem 4.** The following is the DNA sequence of Gene Z that you want to amplify by PCR.

a) Which set of primers from the options below, would you use for the PCR reaction?

Primer Fw

1: 5'-GAGCTCCACTTATA-3'

2: 5'-AGGTGAATATGAAA-3'

3: 5'-GAAAGTATAAGTGG-3'

4: 5'-TTCATATTCACCTC-3'

Primer Rv

1: 5'-ATAGCTAATGCGCG-3'

2: 5'-CCGCGCATTAGCTA-3'

3: 5'-TATCGATTACGCGC-3'

4: 5'-TGGCGAGTAATCGATA-3'

## **Problem 4.** The following is the DNA sequence of Gene Z that you want to amplify by PCR.

a) Which set of primers from the options below, would you use for the PCR reaction?

Primer Fw	Primer Rv	
1: 5'-GAGCTCCACTTATA-3'	1: 5'-ATAGCTAATGCGCG-3'	
2: 5'-AGGTGAATATGAAA-3'	2: 5'-CCGCGCATTAGCTA-3'	
3: 5'-GAAAGTATAAGTGG-3'	3: 5'-TATCGATTACGCGC-3'	
4: 5'-TTCATATTCACCTC-3'	4: 5'-TGGCGAGTAATCGATA-3	

b) What set of primers from the options below, would you use for the PCR reaction? Order them by preference

```
Set 1: 5'-GGTGAATATGAAAG-3' and 5'-TATCGATTACGCGC-3'
Set 2: 5'-AGGTGAATATGAAA-3' and 5'-TATCGATTACGCGC-3'
Set 3: 5'-TCGAGGTGAATATG-3' and 5'-TATCGATTACGCGC-3'
```

Set 4: 5'-CTCGAGGTGAATAT-3' and 5'-TATCGATTACGCGC-3'

**Problem 5.** What will be the result of a PCR amplification using the following set of primers? Choose from the hypotheses given.

5' ACTTCGTTCGCCGGGGCTCGATCGATATTTGGAAT 3'

### 3' TGAAGCAAGCGGCCCCGAGCTAGCTATAAACCTTA 5'

Primer 1: 5'-GTTC-3'

Primer 2: 5'-GCCC-3'

Primer 3: 5'-TATT-3'

Primer 4: 5'-TAGC-3'

Primer 5: 5'- GGAA-3'

Primer 6: 5 - ATTC-3'

	PCR result
Primers 1+2	
Primers 1+3	
Primers 2+3	
Primers 1+4	
Primers 1+2	

### **Amplification hypotheses:**

- 1. no amplified product because the primers face in opposite directions.
- 2. no amplified product because the primers bind to the same strand.
- 3. A product is amplified.
- 4. No product because one of the primers does not hybridize.