


BMC – Bloco Biologia Molecular - 2022

TP4 – PCR “polymerase chain reaction” (part II)

2 most important parameters


Primer design

Annealing temperature




- Sequence complementary to the chain to be replicated
- 3' ends should point to each other
- Fragment to amplify should not exceed 3Kb

PCR – primer design



- Always display seq 5' - 3' and describe the **name** and **characteristics (number of nt; GC%; Tm)** of the primer.


PCR – primer design



Primer p_xptoFw: 5' - ATGCC....-3' (XX nts; YY GC%; Tm)

- Primer forward (or sense): sequence identical to region 5' of the 5'-3' chain
- Always display seq 5' - 3' and describe the **name** and **characteristics (number of nt; GC%; Tm)** of the primer.

PCR – primer design



Primer p_xptoFw: 5' - ATGCC....-3' (XX nts; YY GC%; Tm)

Primer p_xptoRv: 5' - TTGGCG....-3' (XX nts; YY GC%; Tm)

- Reverse primer (or anti-sense): reversed complementary sequence of the 3' chain region 5'-3'
- Always display seq 5' - 3' and describe the **name** and **characteristics (number of nt; GC%; Tm)** of the primer.

- **Sequence size**
- Annealing temperature
- Melting temperature
- GC Content
- GC clamps
- Possible secondary structures
- Repeats of the same nucleotide
- Avoid secondary structures in the template chain
- Avoid regions of frequent homology

PCR – primer design

Primer size

Primers too short


Probability of the sequence to occur:

$$1 \text{ per } 4^8 = 65\,536 \text{ pbs}$$

Human Genome: 3 200 000 000 pbs

$$3\,200\,000\,000 / 65\,536 = 49\,000$$

Primers with 8 nucleotides in length



PCR – primer design

Primer size

Primers too short

Probability of the sequence to occur:

$$1 \text{ per } 4^8 = 65\,536 \text{ pbs}$$

Human Genome: 3 200 000 000 pbs


$$3\,200\,000\,000 / 65\,536 = 49\,000$$

Probability of the sequence to occur:


$$1 \text{ per } 4^{17} = 17\,179\,869\,184 \text{ pbs}$$

$$3\,200\,000\,000 / 17\,179\,869\,184 < 1$$

Primers with 8 nucleotides in length



Primers with 17 nucleotides in length




- Sequence size
- **Annealing temperature**
- **Melting temperature**
- **GC Content**
- GC clamps
- Possible secondary structures
- Repeats of the same nucleotide
- Avoid secondary structures in the template chain
- Avoid regions of frequent homology

PCR – melting temperature

Melting temperature: temperature at which 50% of the chains in solution are denatured

$$T_m = (4 \times (G+C) + 2 \times (A+T)) \text{ } ^\circ\text{C}$$



Primer sequence: 5' **AGACTCAGAGAGAACCC** 3'

4 Gs 5 Cs 7 As 1 T

$$\begin{aligned} T_m &= (4 \times 9) + (2 \times 8) \\ &= 36 + 16 \\ &= 52^\circ\text{C} \end{aligned}$$

$$T^{\circ} \text{ anneal } (T_a) = T_m - 4^{\circ}\text{C}$$

Low T_a :

non-specific products

High T_a :

avoids non-specific products

lowers the amount of product formed

- The T_m of the 2 primers should not differ by more than 3-5 degrees

PCR – annealing temperature

Annealing temperature determines the specificity of amplification.

The hybridization of 2 single chain DNA molecules depends on the temperature.


Nonspecific hybridization

Low enough – increases the efficiency of the PCR reaction

High enough - minimizes mismatches


(a) Annealing temperature is too high

Primers and templates remain dissociated




(b) Annealing temperature is too low

Mismatched hybrid – not all the correct base pairs have formed



(c) Correct annealing temperature

Priming occurs only at the desired target sites



- Sequence size
- Annealing temperature
- Melting temperature
- GC Content
- **GC clamps**
- **Possible secondary structures**
- **Repeats of the same nucleotide**
- Avoid secondary structures in the template chain
- Avoid regions of frequent homology

GC content approx 50-60% - hybridization efficiency of PCR reaction



GC Clamp

The presence of G or C bases in the last 5 bases at the 3'-OH end of the primer (GC Clamp) promotes higher specificity of the bonds at this extremity and increases duplex stability.

PCR – primer design


Possible secondary structures

Presence of secondary inter- or intramolecular structures can decrease efficiency, or even prevent product formation.


Formation of loops (hairpins)

hybridization efficiency
in pcr reaction 


Hairpin loops:



Formation of dimers

hybridization efficiency in
PCR reaction 

Primer-primer dimers:



ΔG closer to zero as possible

Formation of hairpins/dimers must be avoided especially at 3' (prevents polymerase binding)

Repeats of the same nucleotide

Primers with many repeats of a nucleotide can pair unspecifically and should be avoided.

For example AGCGGGGGATGGGG has 5 repetitions of G and 4 more repetitions of this same nucleotide. Maximum number of repetitions = 4 nt.

AGCGGGGGATGGGG

Dinucleotide repeats

Many dinucleotide repeats (for example ATATATAT) should be avoided as they lead to nonspecific pairings, and formation of secondary structures (primer dimers, hairpins). The maximum reference number for dinucleotide repeats should be 2 (2 repetitions of AT = 4 nts: ATAT).

CGTTATATATACGGATCTCG

Primer

CATTTGCAGCATATATATGCCTAGAGCTTTTCTCAGCT

Template DNA

CGTTATATATACGGATCTCG

Correct hybrid

CGTTATATATACGGATCTCG

CGTTATATATACGGATCTCG


CGTTATATATACGGATCTCG

Mismatch hybrids

- Sequence size
- Annealing temperature
- Melting temperature
- GC Content
- GC clamps
- Possible secondary structures
- Repeats of the same nucleotide
- **Avoid secondary structures in the template chain**
- **Avoid regions of frequent homology**

Avoid secondary structures in the template chain

It is important that the primers do not form stable secondary structures with the template chain.



Avoid homology regions

In order to increase the specificity of primers, avoid regions of frequent homology. Primers designed for a specific area should not hybridize other zones at the same locus, or at different loci. Typically, primers are designed and aligned by BLAST to test specificity with the target.

Primer design video:

https://www.youtube.com/watch?v=c-f1H07D_70